

Rabbit Anti-Mouse FcγRI/Cd64 Monoclonal Antibody For Fc - P26151

Item Number: CM0009625

Introduction

Rabbit anti-Mouse FcγRI/CD64 monoclonal antibody validated for flow cytometry (FC). Specific for mouse CD64, a high affinity IgG Fc receptor on macrophages. Conjugated to for direct detection.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse FcγRI/CD64 mAb
Remarks/Alias	CD64; IGGHAFC; FcγgammaRI
Species	Mouse
GenelD (Human)	2209
GenelD	14129
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-297 of mouse FcγRI/CD64 (NP_034316.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P26151
Protein Weight	44kDa
Shipping	Ice bag