

Rabbit Anti-Mouse FcγRI/Cd64 Monoclonal Antibody For Fc - P26151

Item Number: CM0005812

Introduction

Rabbit anti-Mouse FcγRI/CD64 monoclonal antibody for flow cytometry (FC). Detects high-affinity IgG Fc receptor on macrophages, essential for innate and adaptive immunity research. Suitable for mouse immunology studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse FcγRI/CD64 mAb
Remarks/Alias	CD64; IGGHAFC; FcγgammaRI
Species	Mouse
GeneID (Human)	2209
GeneID	14129
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-297 of mouse FcγRI/CD64 (NP_034316.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P26151
Protein Weight	44kDa
Shipping	Ice bag