

Rabbit Anti-Mouse FcγRI/Cd64 Monoclonal Antibody For Flow Cytometry - P26151

Item Number: CM0006132

Introduction

Rabbit monoclonal antibody targeting mouse FcγRI/CD64 (SWISS: P26151), conjugated with . Validated for flow cytometry, cross-reacts with mouse, 44 kDa. Ideal for innate and adaptive immune research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse FcγRI/CD64 mAb
Remarks/Alias	CD64; IGGHAFC; FcγgammaRI
Species	Mouse
GenelD (Human)	2209
GenelD	14129
Immunogen	Recombinant protein - Recombinant fusion protein containing a sequence corresponding to amino acids 25-297 of mouse FcγRI/CD64 (NP_034316.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P26151
Protein Weight	44kDa
Shipping	Ice bag