

Rabbit Anti-Mouse FcγRI/CD64 Monoclonal Antibody For Flow Cytometry - P26151

Item Number: CM0009768



Introduction

High affinity -conjugated rabbit monoclonal antibody against mouse FcγRI/CD64, validated for flow cytometry. Ideal for macrophage identification in innate and adaptive immunity studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse FcγRI/CD64 mAb
Remarks/Alias	CD64; IGGHAFC; FcγgammaRI
Species	Mouse
GenelD (Human)	2209
GenelD	14129
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P26151
Protein Weight	44kDa
Shipping	Ice bag