

Rabbit Anti-Mouse FcγRI/Cd64 Monoclonal Antibody For Flow Cytometry - P26151

Item Number: CM0005241

Introduction

-conjugated rabbit monoclonal antibody targeting mouse FcγRI/CD64, validated for flow cytometry. This high-affinity IgG Fc receptor is macrophage-specific, playing key roles in innate and adaptive immunity. For research use only.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse FcγRI/CD64 mAb
Remarks/Alias	CD64; IGGHAFC; FcγgammaRI
Species	Mouse
Gene ID (Human)	2209
Gene ID	14129
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-297 of mouse FcγRI/CD64 (NP_034316.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Mouse
SWISS	P26151
Protein Weight	44kDa
Shipping	Ice bag