

Rabbit Anti-Mouse Flt3/Cd135 Monoclonal Antibody For Fc - Q00342

Item Number: CM0010647

Introduction

-conjugated rabbit monoclonal antibody targeting mouse Flt3/CD135, validated for flow cytometry. Ideal for hematopoietic stem and progenitor cell and dendritic cell differentiation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Flt3/CD135 mAb
Remarks/Alias	Flk2; Ly72; wmf1; CD135; Flk-2; Flt-3; B230315G04
Species	Mouse
Gene ID (Human)	2322
Gene ID	14255
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00342
Protein Weight	114kDa
Shipping	Ice bag