

Rabbit Anti-Mouse Foxp3 Mab For Fc (Intra) - Q99Jb6

Item Number: CM0011969



Introduction

Rabbit anti-Mouse FOXP3 monoclonal antibody for intracellular flow cytometry. Reacts with mouse FOXP3 (scurfin), a transcriptional regulator essential for Treg development. Validated for mouse samples. - Q99JB6

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse FOXP3 mAb
Remarks/Alias	sf; JM2; scurfin
Species	Mouse
Gene ID (Human)	50943
Gene ID	20371
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse Foxp3 (NP_473380.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q99JB6
Protein Weight	47kDa
Shipping	Ice bag