


Rabbit Anti-Mouse I-A Monoclonal Antibody For Fc - P14483

Item Number: CM0005772



Introduction

-conjugated rabbit monoclonal antibody specific for mouse I-A MHC class II alloantigens (H2-Ab1). Validated for flow cytometry, recognizing the extracellular domain. UniProt: P14483.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse I-A mAb
Remarks/Alias	H-2I; I-Ad, -Ed; I-Eb and I-Ek MHC class II alloantigens; Ia Ag; MHC II
Species	Mouse
GeneID (human)	474382
GeneID	14961
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 30-219 of mouse H2-Ab1 (NP_996988.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P14483
Protein Weight	30kDa
Shipping	Ice bag