

Rabbit Anti-Mouse Icos/Cd278 Monoclonal Antibody For Flow Cytometry - Q9Wvs0

Item Number: CM0006177

Introduction

-conjugated rabbit monoclonal antibody targeting mouse ICOS/CD278, validated for flow cytometry. A costimulatory receptor critical for T-cell responses and antibody class switching.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse ICOS/CD278 mAb
Remarks/Alias	H4; CCLP; AILIM; CRP-1; Ly115
Species	Mouse
GeneID (human)	29851
GeneID	54167
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9WVS0
Protein Weight	22KDa
Shipping	Ice bag