

Rabbit Anti-Mouse Icos/Cd278 Monoclonal Antibody For Flow Cytometry - Q9Wvs0

Item Number: CM0010145

Introduction

-conjugated rabbit monoclonal antibody targeting mouse ICOS/CD278 for flow cytometry. Detects a key T-cell costimulatory receptor involved in immune responses. Cross-reacts with mouse. Protein weight ~22 KDa. Q9WVS0.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse ICOS/CD278 mAb
Remarks/Alias	H4; CCLP; AILIM; CRP-1; Ly115
Species	Mouse
GeneID (Human)	29851
GeneID (Mouse)	54167
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 21-142 of mouse ICOS/CD278 (NP_059508.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9WVS0
Protein Weight	22 KDa
Shipping	Ice bag