

Rabbit Anti-Mouse Ifn-Gamma Monoclonal Antibody For Fc (Intra) - P01580

Item Number: CM0010520

Introduction

-conjugated rabbit monoclonal antibody specific to mouse interferon gamma (IFNG, P01580). Designed for intracellular flow cytometry, cross-reacts with mouse, and supports immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IFN-gamma mAb
Remarks/Alias	Ifg; If2f; IFN-g
Species	Mouse
GeneID (Human)	3458
GeneID	15978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-155 of mouse IFN-gamma (NP_032363.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P01580
Protein Weight	18 kDa
Shipping	Ice bag