

Rabbit Anti-Mouse Ifn-Gamma Monoclonal Antibody For Flow Cytometry (Intracellular) - P01580

Item Number: CM0012349

Introduction

-conjugated rabbit anti-mouse IFN-gamma monoclonal antibody for intracellular flow cytometry. Specific for mouse IFN-gamma, this antibody supports immunological research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IFN-gamma mAb
Remarks/Alias	Ifg; If2f; IFN-g
Species	Mouse
GeneID (Human)	3458
GeneID	15978
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-155 of mouse IFN-gamma (NP_032363.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P01580
Protein Weight	18 kDa
Shipping	Ice bag