

Rabbit Anti-Mouse Iga (Fc) Mab

Item Number: CM0004996

Introduction

-conjugated rabbit monoclonal antibody targeting mouse IgA (Fc region). Designed for intracellular flow cytometry (FC). Validated in mouse samples with recombinant immunogen spanning amino acids 99-344. Ideal for mucosal immunity research.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Mouse IgA (Fc) mAb
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-344 of mouse IgA
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
Shipping	Ice bag