

Rabbit Anti-Mouse Iga (Fc) Monoclonal Antibody

Item Number: CM0004781

Introduction

High-quality -conjugated rabbit monoclonal antibody specifically recognizing mouse IgA Fc region. Ideal for intracellular flow cytometry (FC intra) in immunology and mucosal immunity studies. Reliable IVD raw material from CamelBio.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Mouse IgA (Fc) mAb
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-344 of mouse IgA
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
Shipping	Ice bag