


Rabbit Anti-Mouse IgE Monoclonal Antibody For Fc - P06336

Item Number: CM0005253



Introduction

-conjugated rabbit monoclonal antibody targeting the constant region of mouse IgE. Suitable for intracellular flow cytometry. Reacts specifically with mouse IgE. Molecular weight approximately 47 kDa.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Mouse IgE mAb
Remarks/Alias	IGHE, Ig epsilon chain C region
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 83-421 of mouse IgE (P06336)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P06336
Protein Weight	47kDa
Shipping	Ice bag