


Rabbit Anti-Mouse IgE Monoclonal Antibody For Flow Cytometry (Intracellular) - P06336

Item Number: CM0011739



Introduction

-conjugated rabbit monoclonal antibody targeting mouse IgE heavy chain constant region. Validated for intracellular flow cytometry. Ideal for detecting IgE-expressing B cells and plasma cells in mouse models of allergy and immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IgE mAb
Remarks/Alias	IGHE, Ig epsilon chain C region
Species	Mouse
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 83-421 of mouse IgE (P06336)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P06336
Protein Weight	47kDa
Shipping	Ice bag