


Rabbit Anti-Mouse IgE Monoclonal Antibody For Flow Cytometry (Intracellular) - P06336

Item Number: CM0008960



Introduction

-conjugated rabbit monoclonal antibody against mouse IgE constant region (P06336). Validated for intracellular flow cytometry. Suitable for studying IgE-mediated immunity in allergic and parasitic diseases.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IgE mAb
Remarks/Alias	IGHE, Ig epsilon chain C region
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 83-421 of mouse IgE (P06336).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P06336
Protein Weight	47kDa
Shipping	Ice bag