


Rabbit Anti-Mouse IgG1 Mab For Fc - P01868

Item Number: CM0011425



Introduction

Rabbit anti-Mouse IgG1 monoclonal antibody for flow cytometry. Recombinant protein immunogen, specific for mouse IgG1. Ideal for diagnostic assay development and research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IgG1 mAb
Remarks/Alias	IgG1
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 98-324 of mouse IgG1 (P01868)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01868
Protein Weight	36kDa
Shipping	Ice bag