


Rabbit Anti-Mouse IgG1 Monoclonal Antibody For Flow Cytometry - P01868

Item Number: CM0009721



Introduction

-conjugated rabbit monoclonal antibody targeting mouse IgG1 constant region. Validated for flow cytometry (FC) with mouse samples. Ideal for isotype-specific detection in immunological research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IgG1 mAb
Remarks/Alias	IgG1
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 98-324 of mouse IgG1 (P01868)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01868
Protein Weight	36kDa
Shipping	Ice bag