


Rabbit Anti-Mouse IgG1 Monoclonal Antibody For Flow Cytometry - P01868

Item Number: CM0006003



Introduction

High-affinity -conjugated rabbit anti-mouse IgG1 monoclonal antibody for flow cytometry (FC). Recombinant immunogen targets the constant region, validated for mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IgG1 mAb
Remarks/Alias	IgG1
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 98-324 of mouse IgG1(P01868)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01868
Protein Weight	36kDa
Shipping	Ice bag