

Rabbit Anti-Mouse Igg2B Monoclonal Antibody For Fc - P01867-2

Item Number: CM0008741

Introduction

-conjugated rabbit monoclonal antibody targeting mouse IgG2b (Igh-3). Validated for flow cytometry (FC) with mouse samples. Ideal for immune profiling and isotype detection.

[Learn More](#)

Parameter	Value
Product Name	Rabbit Anti-Mouse IgG2b mAb
Remarks/Alias	IgG2b, Igh-3, gamma2b
Species	Mouse
GeneID	16016
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 14-113 of mouse IgG2b (P01867-2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01867-2
Protein Weight	37kDa
Shipping	Ice bag