

Rabbit Anti-Mouse IgM Mab - P01872

Item Number: CM0010872



Introduction

-conjugated rabbit monoclonal antibody targeting the constant region of mouse IgM heavy chain. Validated for flow cytometry (FC) detection of mouse IgM-expressing B cells, supporting immunoglobulin and humoral immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IgM mAb
Remarks/Alias	Igm; muH; Igh6; Igh-6; Igh-M
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 105-454 of mouse IgM (P01872)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01872
Protein Weight	50kDa
Shipping	Ice bag