

Rabbit Anti-Mouse IgM Monoclonal Antibody For Fc - P01872

Item Number: CM0009835



Introduction

-conjugated rabbit anti-Mouse IgM monoclonal antibody validated for flow cytometry (FC) in mouse. Recombinant immunogen covers amino acids 105–454 of the IgM heavy chain constant region (P01872). Ideal for B cell immunophenotyping and humoral immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IgM mAb
Remarks/Alias	Igm; muH; Igh6; Igh-6; Igh-M
Species	Mouse
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 105-454 of mouse IgM (P01872)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01872
Protein Weight	50kDa/53kDa
Shipping	Ice bag