

Rabbit Anti-Mouse IgM Monoclonal Antibody For Flow Cytometry - P01872

Item Number: CM0004874



Introduction

-conjugated rabbit monoclonal antibody targeting mouse IgM heavy constant region (P01872), validated for flow cytometry. Ideal for B cell immunophenotyping, immune response studies, and mucosal immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IgM mAb
Remarks/Alias	Igm; muH; Igh6; Igh-6; Igh-M
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 105-454 of mouse IgM (P01872)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01872
Protein Weight	50kDa
Shipping	Ice bag