


Rabbit Anti-Mouse IL-12 P40 Monoclonal Antibody For Flow Cytometry (Intracellular) - P43432

Item Number: CM0010852



Introduction

-conjugated rabbit monoclonal antibody against mouse IL-12 p40, validated for intracellular flow cytometry. Recognizes the p40 subunit common to IL-12 and IL-23, key regulators of immunity and inflammation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL-12 p40 mAb
Remarks/Alias	p40; IL-12b; IL12p40; IL-12p40
Species	Mouse
GeneID (Human)	3593
GeneID	16160
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-335 of mouse IL-12 p40(NP_001290173.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P43432
Protein Weight	38kDa
Shipping	Ice bag