

Rabbit Anti-Mouse IL-13 Monoclonal Antibody For Flow Cytometry (Intracellular) - P20109

Item Number: CM0005602

Introduction

-conjugated rabbit monoclonal antibody targeting mouse Interleukin-13 (IL-13). Designed for intracellular flow cytometry, suitable for studies on allergic inflammation, cytokine signaling, and immune response.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL-13 mAb
Remarks/Alias	IL-13
Species	Mouse
GenelD (human)	3596
GenelD	16163
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P20109
Protein Weight	14kDa
Shipping	Ice bag