

Rabbit Anti-Mouse IL-13 Monoclonal Antibody For Intracellular Flow Cytometry - P20109

Item Number: CM0004134

Introduction

-conjugated rabbit monoclonal antibody targeting mouse IL-13, validated for intracellular flow cytometry. Reacts with mouse, 14 kDa cytokine involved in allergic inflammation and immune response. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL-13 mAb
Remarks/Alias	IL-13
Species	Mouse
GeneID (human)	3596
GeneID	16163
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-131 of mouse IL-13 (NP_032381.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P20109
Protein Weight	14kDa
Shipping	Ice bag