

Rabbit Anti-Mouse IL-17A Monoclonal Antibody For Flow Cytometry - Q62386

Item Number: CM0011344

Introduction

-conjugated rabbit monoclonal antibody against mouse IL-17A, validated for intracellular flow cytometry. Recognizes the 18 kDa cytokine, suitable for Th17 and innate immunity research in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL-17A mAb
Remarks/Alias	IL17; Ctlα8; IL-17; Ctlα-8; IL-17A
Species	Mouse
GeneID (Human)	3605
GeneID	16171
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-158 of mouse IL17A. (NP_034682.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q62386
Protein Weight	18 kDa
Shipping	Ice bag