

Rabbit Anti-Mouse IL-9 Mab - P15247

Item Number: CM0004923



Introduction

-conjugated rabbit monoclonal antibody against mouse interleukin-9 (IL-9, P40) for intracellular flow cytometry. Validated for mouse samples, targets a key Th2 cytokine involved in immune response and mast cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL-9 mAb
Remarks/Alias	P40; IL-9
Species	Mouse
Gene ID (Human)	3578
Gene ID	16198
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-114 of mouse IL-9 (NP_032399.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P15247
Protein Weight	16kDa
Shipping	Ice bag