

Rabbit Anti-Mouse IL1R1/Cd121A Mab - P13504

Item Number: CM0008246

Introduction

-conjugated rabbit monoclonal antibody against mouse IL1R1/CD121a. Validated for flow cytometry (FC). Reacts specifically with mouse IL-1 receptor type 1. Ideal for immunological and neuroinflammatory studies.

[Learn More](#)

Specification	Details
Product Name	Rabbit anti-Mouse IL1R1/CD121a mAb
Remarks/Alias	IL-1R; CD121a; CD121b; IL-1R1; IL1r-1; IL-1R-1; IL-1RT1; IL-1RT-1; IL-1R-alpha
Species	Mouse
GeneID (human)	3554
GeneID	16177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-338 of mouse IL1R1/CD121a (NP_032388.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13504
Protein Weight	67kDa
Shipping	Ice bag