

Rabbit Anti-Mouse IL1R1/Cd121A Monoclonal Antibody For Flow Cytometry - P13504

Item Number: CM0010088

Introduction

-conjugated rabbit monoclonal antibody against mouse IL-1 receptor type 1 (IL1R1/CD121a) for flow cytometry. Validated in mouse samples. Ideal for IL-1 signaling and immune research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL1R1/CD121a mAb
Remarks/Alias	IL-1R; CD121a; CD121b; IL-1R1; IL1r-1; IL-1R-1; IL-1RT1; IL-1RT-1; IL-1R-alpha
Species	Mouse
GeneID (Human)	3554
GeneID	16177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-338 of mouse IL1R1/CD121a (NP_032388.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13504
Protein Weight	67kDa
Shipping	Ice bag