

Rabbit Anti-Mouse IL22 Monoclonal Antibody For Intracellular Flow Cytometry - Q9Jjy9

Item Number: CM0009621



Introduction

-conjugated rabbit monoclonal antibody specific to mouse IL-22, validated for intracellular flow cytometry (FC intra). Ideal for research on cytokine signaling, tissue regeneration, and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL22 mAb
Remarks/Alias	IL-22; If2b1; Itif; IL-22a; ILT1Fa
Species	Mouse
Gene ID (Human)	50616
Gene ID (Mouse)	50929
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q9Jjy9
Protein Weight	20kDa
Shipping	Ice bag