

Rabbit Anti-Mouse IL22 Monoclonal Antibody For Intracellular Flow Cytometry - Q9Jjy9

Item Number: CM0009506



Introduction

-conjugated rabbit monoclonal antibody against mouse IL22, validated for intracellular flow cytometry. Ideal for investigating IL-22-mediated epithelial regeneration, barrier function, and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL22 mAb
Remarks/Alias	IL-22; If2b1; Itif; IL-22a; ILT1Fa
Species	Mouse
GeneID (Human)	50616
GeneID	50929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-179 of mouse IL22 (NP_058667.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q9JJY9
Protein Weight	20kDa
Shipping	Ice bag