

Rabbit Anti-Mouse Integrin Alpha 1/Cd49A Monoclonal Antibody For Fc - Q3V3R4

Item Number: CM0011060

Introduction

-conjugated rabbit anti-mouse Integrin alpha 1/CD49a monoclonal antibody for flow cytometry. Recognizes laminin/collagen receptor involved in cell adhesion and growth regulation. Reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Integrin alpha 1/CD49a mAb
Remarks/Alias	Vla1; CD49A; E130012M19Rik
Species	Mouse
Gene ID (Human)	3672
Gene ID	109700
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-1141 of mouse Integrin alpha 1/CD49a (NP_001028400.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q3V3R4
Protein Weight	130kDa
Shipping	Ice bag