

Rabbit Anti-Mouse Integrin Alpha 2/Cd49B Monoclonal Antibody For Flow Cytometry - Q62469

Item Number: CM0012642

Introduction

-conjugated rabbit monoclonal antibody against mouse Integrin alpha 2/CD49b (DX5, GPIa, CD49B). Validated for flow cytometry (FC). Integrin alpha-2/beta-1 is a collagen receptor involved in cell adhesion and extracellular matrix organization.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Integrin alpha 2/CD49b mAb
Remarks/Alias	DX5; GPIa; CD49B
Species	Mouse
GeneID (Human)	3673
GeneID (Mouse)	16398
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 168-366 of mouse Integrin alpha 2/CD49b (NP_032422.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q62469
Protein Weight	129kDa
Shipping	Ice bag