

Rabbit Anti-Mouse Integrin Alpha 2/Cd49B Monoclonal Antibody For Flow Cytometry - Q62469

Item Number: CM0009174



Introduction

Rabbit monoclonal antibody specifically binding mouse Integrin alpha-2/CD49b, validated for flow cytometry. Recognizes a collagen receptor crucial for cell adhesion and platelet function. Suitable for mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Integrin alpha 2/CD49b mAb
Remarks/Alias	DX5; GPIa; CD49B
Species	Mouse
GenelD (human)	3673
GenelD	16398
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 168-366 of mouse Integrin alpha 2/CD49b (NP_032422.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q62469
Protein Weight	129kDa
Shipping	Ice bag