

Rabbit Anti-Mouse Jam-1/Cd321 Monoclonal Antibody For Fc - O88792

Item Number: CM0004681



Introduction

-conjugated rabbit monoclonal antibody targeting mouse JAM-1/CD321 (UniProt O88792). Designed for flow cytometry, this antibody aids research on tight junctions, leukocyte transmigration, and platelet activation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse JAM-1/CD321 mAb
Remarks/Alias	JAM; Jcam; JAM-1; JAM-A; Jcam1; Ly106; ESTM33; 9130004G24; JAM-A/CD321
Species	Mouse
GeneID (Human)	50848
GeneID	16456
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-242 of mouse JAM-1/CD321 (NP_766235.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O88792
Protein Weight	32kDa
Shipping	Ice bag