

Rabbit Anti-Mouse Jam-1/Cd321 Monoclonal Antibody For Flow Cytometry - O88792

Item Number: CM0008138



Introduction

-conjugated rabbit monoclonal antibody against mouse JAM-1/CD321, validated for flow cytometry (FC). Recognizes the 32kDa junctional adhesion molecule A involved in tight junction formation and leukocyte transmigration.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Mouse JAM-1/CD321 mAb
Remarks/Alias	JAM; Jcam; JAM-1; JAM-A; Jcam1; Ly106; ESTM33; 9130004G24; JAM-A/CD321
Species	Mouse
GeneID (Human)	50848
GeneID	16456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-242 of mouse JAM-1/CD321 (NP_766235.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O88792
Protein Weight	32kDa
Shipping	Ice bag