


Rabbit Anti-Mouse Lamp1/Cd107A Monoclonal Antibody For Flow Cytometry - P11438

Item Number: CM0004936



Introduction

-conjugated rabbit anti-mouse LAMP1/CD107a monoclonal antibody for flow cytometry. Targets a lysosomal membrane glycoprotein involved in autophagy, NK cell cytotoxicity, and metastasis. Recombinant immunogen corresponding to amino acids 208-363.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse LAMP1/CD107a mAb
Remarks/Alias	P2B; Perk; LGP-A; CD107a; Lamp-1; LGP-120
Species	Mouse
GeneID (Human)	3916
GeneID (Mouse)	16783
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 208-363 of mouse LAMP1/CD107a (NP_034814.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P11438
Protein Weight	43kDa
Shipping	Ice bag