


Rabbit Anti-Mouse Lamp1/Cd107A Monoclonal Antibody For If/Icc, Fc - P11438

Item Number: CM0011296



Introduction

-conjugated rabbit monoclonal antibody targeting mouse LAMP1/CD107a for immunofluorescence and flow cytometry. Ideal for lysosome biogenesis, autophagy, and NK cell cytotoxicity studies. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse LAMP1/CD107a mAb
Remarks/Alias	P2B; Perk; LGP-A; CD107a; Lamp-1; LGP-120
Species	Mouse
GeneID (Human)	3916
GeneID (Mouse)	16783
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 208-363 of mouse LAMP1/CD107a (NP_034814.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	P11438
Protein Weight	43kDa
Shipping	Ice bag