

Rabbit Anti-Mouse Ly-49D Monoclonal Antibody For Fc - Q60651

Item Number: CM0011832

	Introduction -conjugated rabbit monoclonal antibody targeting mouse Ly-49D (KLRA4), validated for flow cytometry analysis of NK cell receptor expression. Reacts with mouse. Ideal for immunology research.
--	---

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Ly-49D mAb
Remarks/Alias	Chok; Ly49d; Klra32; Klra33; Ly-49d; ly49r<129>
Species	Mouse
Gene ID	16635
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 66-263 of mouse Ly-49D (NP_001239506.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q60651
Protein Weight	31kDa
Shipping	Ice bag