

Rabbit Anti-Mouse Ly-49D Monoclonal Antibody For Flow Cytometry - Q60651

Item Number: CM0011613



Introduction

-conjugated rabbit monoclonal antibody targeting mouse Ly-49D (Killer cell lectin-like receptor 4). Validated for flow cytometry, ideal for NK cell receptor and immune profiling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Ly-49D mAb
Remarks/Alias	Chok; Ly49d; Klra32; Klra33; Ly-49d; ly49r<129>
Species	Mouse
GenelD	16635
Immunogen	Recombinant fusion protein containing amino acids 66-263 of mouse Ly-49D (NP_001239506.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Mouse
SWISS	Q60651
Protein Weight	31kDa
Shipping	Ice bag