

Rabbit Anti-Mouse Ly-6A/E (Sca-1) Mab For Fc - P08648

Item Number: CM0010275



Introduction

-conjugated rabbit monoclonal antibody targeting mouse Ly-6A/E (Sca-1) for flow cytometry. Validated for mouse samples, ideal for stem cell and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Ly-6A/E (Sca-1) mAb
Remarks/Alias	TAP; Sca1; Sca-1; Ly-6A.2; Ly-6A/E; Ly-6E.1
Species	Mouse
GeneID	110454
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P08648
Protein Weight	114kDa
Shipping	Ice bag