

Rabbit Anti-Mouse Ly-6A/E (Sca-1) Monoclonal Antibody For Flow Cytometry - P08648

Item Number: CM0009946



Introduction

PerCP/Cyanine5.5-conjugated rabbit monoclonal antibody targeting mouse Ly-6A/E (Sca-1). Validated for flow cytometry, reacts with mouse, 115kDa. Suitable for stem cell and hematopoietic research.

[Learn More](#)

Parameter	Value
Product Name	PerCP/Cyanine5.5 Rabbit anti-Mouse Ly-6A/E (Sca-1) mAb
Remarks/Alias	TAP; Sca1; Sca-1; Ly-6A.2; Ly-6A/E; Ly-6E.1
Species	Mouse
Gene ID	110454
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-111 of mouse Ly-6A/E (Sca-1) (NP_034868.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P08648
Protein Weight	115kDa
Shipping	Ice bag