

Rabbit Anti-Mouse Pd-L2/Cd273/B7-Dc Monoclonal Antibody For Flow Cytometry - Q9Wul5

Item Number: CM0004794



Introduction

-conjugated rabbit anti-mouse PD-L2 (CD273/B7-DC) monoclonal antibody for flow cytometry. Validate PD-L2 expression in immune tolerance and cancer research. Specific for mouse samples. - Q9WUL5

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse PD-L2/CD273/B7-DC mAb
Remarks/Alias	Btdc; B7-DC; PD-L2; F730015O22Rik
Species	Mouse
Gene ID (Human)	80380
Gene ID	58205
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-123 of mouse PD-L2/CD273/B7-DC (NP_067371.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9WUL5
Protein Weight	28kDa
Shipping	Ice bag