

## Rabbit Anti-Mouse Rank/Cd265 Monoclonal Antibody For Fc - O35305

Item Number: CM0008760



### Introduction

-conjugated rabbit monoclonal antibody against mouse RANK/CD265 (TNFRSF11A) for flow cytometry (FC). Recognizes mouse RANK, essential for osteoclastogenesis and immune regulation. Validated for mouse samples.

[Learn More](#)

| Parameter        | Value                                                                                                                                       |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse RANK/CD265 mAb                                                                                                            |
| Remarks/Alias    | OFE; ODFR; Rank; Ly109; mRANK; TRANCE-R                                                                                                     |
| Species          | Mouse                                                                                                                                       |
| GeneID (human)   | 8792                                                                                                                                        |
| GeneID           | 21934                                                                                                                                       |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 31-213 of mouse RANK/CD265 (NP_033425.3). |
| Source           | Rabbit                                                                                                                                      |
| Category         | Monoclonal Antibodies                                                                                                                       |
| Application      | FC                                                                                                                                          |
| Cross Reactivity | Mouse                                                                                                                                       |
| SWISS            | O35305                                                                                                                                      |
| Protein Weight   | 67kDa                                                                                                                                       |
| Shipping         | Ice bag                                                                                                                                     |