

Rabbit Anti-Mouse Rank/Cd265 Monoclonal Antibody For Flow Cytometry - O35305

Item Number: CM0004422



Introduction

-conjugated rabbit monoclonal antibody against mouse RANK/CD265 (TNFRSF11A), validated for flow cytometry. Suitable for studying osteoclast differentiation and T-cell-dendritic cell interactions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse RANK/CD265 mAb
Remarks/Alias	OFE; ODFR; Rank; Ly109; mRANK; TRANCE-R
Species	Mouse
GeneID (Human)	8792
GeneID	21934
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-213 of mouse RANK/CD265 (NP_033425.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O35305
Protein Weight	67kDa
Shipping	Ice bag