


Rabbit Anti-Mouse S100A9 Monoclonal Antibody For Flow Cytometry (Intracellular) - P31725

Item Number: CM0008501



Introduction

-conjugated rabbit monoclonal antibody targeting mouse S100A9 (Calgranulin-B/MRP14) for intracellular flow cytometry. Validated in mouse. S100A9 functions in inflammatory response and immune regulation. 13 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse S100A9 mAb
Remarks/Alias	p14; Cagb; GAGB; L1Ag; BEE22; MRP14; 60B8Ag
Species	Mouse
GeneID (Human)	6280
GeneID	20202
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of mouse S100a9 (NP_033140.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P31725
Protein Weight	13kDa
Shipping	Ice bag