

Rabbit Anti-Mouse T-Bet/Tbx21 Monoclonal Antibody For Intracellular Flow Cytometry - Q9Jkd8

Item Number: CM0008968



Introduction

-conjugated rabbit monoclonal antibody targeting mouse T-bet/Tbx21, a Th1 lineage transcription factor. Validated for intracellular flow cytometry, it enables specific T-bet detection in mouse T cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse T-bet/Tbx21 mAb
Remarks/Alias	TBT1; Tbet; Tblym
Species	Mouse
GeneID (human)	30009
GeneID	57765
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 321-530 of mouse T-bet/Tbx21 (NP_062380.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q9JKD8
Protein Weight	58 kDa
Shipping	Ice bag