

Rabbit Anti-Mouse Vcam-1/Cd106 Monoclonal Antibody For Fc - P29533

Item Number: CM0006129

Introduction

Rabbit anti-Mouse VCAM-1/CD106 monoclonal antibody validated for flow cytometry (FC) in mouse samples. Enables investigation of leukocyte adhesion and vascular inflammation. UniProt: P29533.

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Mouse VCAM-1/CD106 mAb
Remarks/Alias	CD106; Vcam-1
Species	Mouse
GenelD (human)	7412
GenelD	22329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-310 of mouse VCAM-1/CD106 (NP_035823.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P29533
Protein Weight	38kDa/81kDa
Shipping	Ice bag